

Features

- 55" 4K Inwall console
- Integrates with CREA Essential as secondary display, controlled via CREA Essential
- 4 x 4K inputs with multiview capabilities - PIP PAP QUAD
- IP 65 rated Metal powder coated inwall enclosure with AV systems flush mounted inside OT wall
- Choice of medical grade and commercial grade displays

Specifications

Display Size	55"		
I/O	Inputs	USB 2.0 Type A	1
		HDMI 2.0	2
		LAN	1
	Output	HDMI 2.0	1 (optional)
Physical Characteristics	Dimensions (W x H x D)mm		1445 x 1015 x 153
	Weight		70 -100 kgs (based on display option)
Mounting	In wall		

Option 1	55" Medical Grade display panel	Diagonal size	55"
		Resolution	3840 x 2160 (16:9)
		Brightness	350 nit
		Contrast ratio	1,000:1
		Viewing angle	178 / 178
		DICOM Part 14 Compliant	
		True flat design	
		Displayable colors	1.07 Billion (10-bit)

Option 3	55" Commercial Grade display panel	Diagonal size	55"
		Resolution	3840 x 2160 (16:9)
		Brightness	350 nit
		Contrast ratio	4,000:1
		Displayable colors	1.07 Billion (10-bit)

Power	AC 100 - 240 V, Max 250 W	
Environment	Operating Temperature	0 ~ 40°C / 32 ~ 122° F